

Class B Sterilisers

Meet the NEW Midmark Newmed Class B Steriliser

With fast cycles suitable for handpieces

Good looking and elegant on the outside, but reliable and uncomplicated on the inside, the Midmark Newmed Class B Steriliser is your perfect partner for instrument processing.

Our Italian design provides sleek, stylish lines, to enhance the look of even the parts of your surgery unvisited by patients, whilst maintaining easily wipeable surfaces for decontamination and integral data capture for effortless record keeping.

Our experienced, forward thinking team of engineers have taken the opportunity to design a steriliser range with a faster printer, faster door closure mechanism, improved data recording functions, a more intuitive user interface and more intelligent diagnostics to simplify maintenance procedures and maximise uptime.

Whether you use your steriliser all day every day or less frequently, the Midmark Newmed is the reliable, cost effective unit for your practice, processing unwrapped instruments for immediate use or pouched items which can be stored until required.*

*maximum storage time depends on local regulations and best practices



Midmark EMEA Ltd. Beech House, First Floor, East Wing, Ancells Business Park, Fleet, Hampshire, GU51 2UN, UK.
Tel: +44 1252 360940 Fax: +44 1252 360941

www.midmark.com



Class B Sterilisers

Designed for perfect sterilisation of your materials, thanks to our precise combination of vacuum, temperature and drying parameters


Table of Cycles

Materials and instruments	Program	Vacuum phases	Total time*
Delicate hollow** and rustproof hollow (wrapped)	121* Hollow wrapped	4	28 minutes
Rustproof hollow (wrapped)	134* Hollow wrapped	4	15 minutes
Solid made of rubber and delicate solids (wrapped)	121* Hollow wrapped	2	28 minutes
Rustproof solid (wrapped)	134* Hollow wrapped	2	15 minutes
Rustproof solid and hollow (wrapped)	134* Prion	4	30 minutes
Delicate, hollow and porous (wrapped)	121* Porous	4	33 minutes
Solid, rustproof, hollow, and porous (wrapped)	134* Porous	4	20 minutes
Solid made of rubber and delicate (unwrapped)	121* Rapid	2	22 minutes
Solid rubber (unwrapped)	134* Rapid	2	9 minutes
Hollow made of rubber and stainless steel (unwrapped)	134* Open hollow	4	9 minutes
Text Cycle for Helix and Bowie & Dick	134* Helix/B&D Test	4	7.5 minutes
Test Cycle for vacuum	<40* Vacuum Test	1	15 minutes

*Includes drying phase, excluding preheating. The preheating may vary from 20-25 minutes depending on the load

* Type A and B hollow instruments, in accordance with EN 13060

*** For Midmark Newmed B23 the total time must be increased by 15 minutes



Midmark EMEA Ltd. Beech House, First Floor, East Wing, Ancells Business Park, Fleet, Hampshire, GU51 2UN, UK.
Tel: +44 1252 360940 Fax: +44 1252 360941

www.midmark.com

